


# Main window



File: Reset

Save path

Exit

Settings: Trigger--

Set timer

Number of pbbes (4-4)

Calibration

Help: View documentation

Sampling frequency  
view trigger - calibration setting threshold  
[Choose calibration file]

# On startup

1.0

No configuration file found. Do you wish to configure it?

Yes

Not now

1.1

Warning!  
You won't be able to run experiments until calibration is done

OK

Cancel

main window

2.0

What kind of triggers will you be using? (Can be changed later)

- ☐ No trigger
- ☐ Mouse sensor
- ☐ Network

OK

Cancel

2.1

Are you sure?

OK

Reset

Cancel

main window

3.0

Make sure 4 mice are connected to the computer

OK

Cancel

4.0

Identify sensors by pressing the button below and move sensors within 10s

1A sensor

Next

Cancel

4.1

Move sensors...

Cancel

